

UK NEQAS EQA

Instructions for processing the UK NEQAS FB, RE and DL EQA Samples on the Abbott Alinity Hq Analysers

Analysis Instructions We would recommend that all UK NEQAS EQA samples for the FBC, Retic and Automated Differential programmes are processed in the normal patient processing mode.

Additional Notes: When reporting the neutrophil concentration and percentage, the IG (Immature Neutrophils) must be included (combined) into the neutrophil concentration.

Contact Information Ludi van Dun, Manager customer engagement and scientific liaison
Lpjm.vandun@abbott.com
+31.653386997

Authorised by



Abbott Diagnostics
Wegalaan 9
2132 JD Hoofddorp
Netherlands

Ludi van Dun
Abbott Diagnostics EMEA