



Ms Yvonne John
Haematology Service Manager
UK NEQAS Haematology
PO Box 14
Watford
WD180FJ
UK

Abbott Diagnostics
Abbott House
Vanwall Business Park
Vanwall Road
Maidenhead
SL6 4XF

10th September 2018

RE: Processing EQA Samples on CELL-DYN 3700

Dear Yvonne,

Thank you for bringing the matter of EQA samples on CELL-DYN 3700 to our attention. We are pleased to be able to offer the following guidance for all CELL-DYN 3700 customers running UK NEQAS samples.

We would recommend that all UK NEQAS FBC EQA samples are processed in the WIC mode, and those results reported.

I hope this provides the guidance you require.

Yours sincerely

PP

Dr. Zainab Mukhtar PhD
Abbott Haematology Medical Affairs