

Mr Paul McTaggart
Data Manager
UK NEQAS Haematology
PO Box 14
Watford
WD18 0FJ
UK

Abbott Diagnostics
Abbott House
Vanwall Business Park
Vanwall Road
Maidenhead
SL6 4XF

16th November 2016

RE: Processing EQA Samples on CELL-DYN Ruby

Dear Paul,

Thank you for bringing the matter of EQA samples on CELL-DYN Ruby to our attentions. We are pleased to be able to offer the following guidance for all CELL-DYN Ruby customers running UK NEQAS samples.

We would recommend that all UK NEQAS EQA samples are processed in the NOC mode, and those results reported.

I hope this provides the guidance you require.

Yours sincerely


pp

Dr. Hans J. Hoffmann PhD

Abbott Haematology Global Scientific Affairs